



Anti-REXO1 polyclonal antibody (DPAB-DC2584)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	REXO1 (REX1, RNA exonuclease 1 homolog (<i>S. cerevisiae</i>)) is a protein-coding gene. GO annotations related to this gene include exonuclease activity and nucleic acid binding. An important paralog of this gene is REXO1L11P.
Immunogen	A synthetic peptide corresponding to human REXO1. The sequence is C-KVREDAKTKR
Source/Host	Goat
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Concentration	0.5 mg/mL
Size	100 µg
Buffer	In Tris saline, pH 7.3 (0.5% BSA, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	REXO1 REX1, RNA exonuclease 1 homolog (S. cerevisiae) [Homo sapiens (human)]
Official Symbol	REXO1
Synonyms	REXO1; REX1, RNA exonuclease 1 homolog (S. cerevisiae); REX1; ELOABP1; EloA-BP1; TCEB3BP1; RNA exonuclease 1 homolog; elongin A binding protein 1; elongin A-binding protein 1; elongin-A-binding protein 1; transcription elongation factor B polypeptide 3 binding protein 1; transcription elongation factor B polypeptide 3-binding protein 1;
Entrez Gene ID	57455
Protein Refseq	NP_065746
UniProt ID	Q8N1G1
Chromosome Location	19p13.3
Pathway	Ribosome biogenesis in eukaryotes;
Function	exonuclease activity; nucleic acid binding;