



Anti-RNF150 (aa 67-176) polyclonal antibody (DPAB-DC2589)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNF150 (ring finger protein 150) is a protein-coding gene. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is ZNRF3.
Immunogen	RNF150 (NP_065775, 67 a.a. ~ 176 a.a) partial recombinant protein with GST tag. The sequence is LHTEKTECGRYGEHSPKQDARGEVVMASAHDRACDPNTKFAAPTRGKNWIALIPKGNC TYRDKIRNAFLQNASAVVIFNVGSNTNETITMPHAGVEDIVAIMIPEPKG
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RNF150 ring finger protein 150 [Homo sapiens (human)]
Official Symbol	RNF150

Synonyms	RNF150; ring finger protein 150; RING finger protein 150;
Entrez Gene ID	57484
Protein Refseq	NP_065775
UniProt ID	Q9ULK6
Chromosome Location	4q31.21
Function	zinc ion binding;