



Anti-BAHCC1 polyclonal antibody (DPAB-DC2599)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	BAHCC1 (BAH domain and coiled-coil containing 1) is a protein-coding gene. GO annotations related to this gene include DNA binding. An important paralog of this gene is TNRC18.
Immunogen	Recombinant protein corresponding to amino acids of recombinant BAHCC1. The sequence is SLGLLCAELRGSGGEPAAKKRSKLSVYAGLQTASVEKAQCKKSSCQGGLAPSVHRVA QLKPKVKSKGLPTGLSSFQQKEATPGGRIREKLSRAK
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	BAHCC1 BAH domain and coiled-coil containing 1 [Homo sapiens (human)]
Official Symbol	BAHCC1
Synonyms	BAHCC1; BAH domain and coiled-coil containing 1; BAHD2; BAH and coiled-coil domain-containing protein 1; BAH domain-containing protein 2; bromo adjacent homology domain-containing protein 2;
Entrez Gene ID	57597
Protein Refseq	NP_001278253
UniProt ID	A0A075B747
Chromosome Location	17q25.3
Function	chromatin binding;