



Anti-PHRF1 (aa 92-189) polyclonal antibody (DPAB-DC2605)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PHRF1 (PHD and ring finger domains 1) is a protein-coding gene. Diseases associated with PHRF1 include systemic lupus erythematosus, and lupus erythematosus. GO annotations related to this gene include RNA polymerase binding and protein domain specific binding. An important paralog of this gene is SCAF1.
Immunogen	KIAA1542 (NP_065952, 92 a.a. ~ 189 a.a) partial recombinant protein with GST tag. The sequence is KLEAAGSFNSDDDAESCPICLNAFRDQAVGTPENCAHYFCLDCIVEWSKNANSCPVDRTL FKCICIRAQFGGKILKKIPVENTKASEEEEDPTFCEVC
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [PHRF1 PHD and ring finger domains 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	PHRF1
Synonyms	PHRF1; PHD and ring finger domains 1; RNF221; PPP1R125; PHD and RING finger domain-containing protein 1; CTD-binding SR-like protein rA9; protein phosphatase 1, regulatory subunit 125;
Entrez Gene ID	57661
Protein Refseq	NP_001273510
UniProt ID	A0A024RCA1
Chromosome Location	11p15.5
Function	RNA polymerase binding; protein domain specific binding; zinc ion binding;