



Anti-C6ORF47 polyclonal antibody (DPAB-DC2616)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C6orf47 (chromosome 6 open reading frame 47) is a protein-coding gene. Diseases associated with C6orf47 include human t-cell leukemia virus type 1, and gastroenteritis.
Immunogen	Recombinant protein corresponding to amino acids of recombinant C6orf47. The sequence is EPSKEEPQVEQLGSKRMDSLKWDQPISSTQESGRLEAGGASPKLRWDHVDSSGGTRRPGVS PEGGLSVPGPGAPLEK
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC-P, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	C6orf47 chromosome 6 open reading frame 47 [Homo sapiens (human)]
Official Symbol	C6ORF47
Synonyms	C6ORF47; chromosome 6 open reading frame 47; G4; NG34; D6S53E; uncharacterized protein C6orf47; protein G4;
Entrez Gene ID	57827
Protein Refseq	NP_067007
UniProt ID	O95873
Chromosome Location	6p21.3