



Anti-ERG28 (N-terminal) polyclonal antibody (DPAB-DC2646)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of mouse ERG28.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	0610007P14Rik RIKEN cDNA 0610007P14 gene [Mus musculus (house mouse)]
Official Symbol	0610007P14RIK
Synonyms	0610007P14RIK; RIKEN cDNA 0610007P14 gene; ORF11; C77855; AU019315; 1190004E09Rik; probable ergosterol biosynthetic protein 28; c14orf1-like protein;

Entrez Gene ID	58520
Protein Refseq	NP_067421
UniProt ID	Q9ERY9
Chromosome Location	12 D2; 12