



Anti-ZFP862-PS polyclonal antibody (DPAB-DC2655)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to recombinant protein GX0784. This antibody detects mZNF862 protein.
Immunogen	Recombinant GST fusion protein corresponding to 189 amino acids of mouse Zfp862.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Zfp862-ps zinc finger protein 862, pseudogene [Mus musculus (house mouse)]
Official Symbol	ZFP862-PS

Synonyms	ZFP862-PS; zinc finger protein 862, pseudogene; Zfp862; AI448302; mKIAA0543; 4732460K03Rik;
Entrez Gene ID	58894
Chromosome Location	6 B2.3; 6