



Anti-HAPLN2 (aa 31-142) polyclonal antibody (DPAB-DC2708)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	HAPLN2 (hyaluronan and proteoglycan link protein 2) is a protein-coding gene. Diseases associated with HAPLN2 include rabies, and malignant glioma, and among its related super-pathways are FAK1 Signaling and MAPK Signaling. GO annotations related to this gene include hyaluronic acid binding. An important paralog of this gene is HAPLN4.
Immunogen	HAPLN2 (NP_068589, 31 a.a. ~ 142 a.a) partial recombinant protein with GST tag. The sequence is HPGPHYLLPPIHEVIHSHRGATATLPCVLGTTTPSYKVRWSKVEPGELRETLILITNGLH ARGYGPLGGRARMRRGHRDASLVIAGVRLEDEGRYRCELINGIEDESVALT
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [HAPLN2 hyaluronan and proteoglycan link protein 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	HAPLN2
Synonyms	HAPLN2; hyaluronan and proteoglycan link protein 2; BRAL1; brain link protein 1; brain link protein-1;
Entrez Gene ID	60484
Protein Refseq	NP_068589
UniProt ID	Q9GZV7
Chromosome Location	1q23.1
Function	hyaluronic acid binding;