



Anti-Maltose Phosphorylase (native) polyclonal antibody (DPAB-DC2757)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified Escherichia coli Maltose Phosphorylase.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Bacteria
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IP, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	LAF_0025 maltose phosphorylase [Lactobacillus fermentum IFO 3956]
Official Symbol	LAF_0025

Synonyms	LAF_0025; maltose phosphorylase;
Entrez Gene ID	6233598
UniProt ID	B2GEX2
Chromosome Location	NC_010610.1 (27477..29723)
Pathway	Metabolic pathways; Starch and sucrose metabolism;