



Anti-PLIN3 (aa 1-100) polyclonal antibody (DPAB-DC2951)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PLIN3 (perilipin 3) is a protein-coding gene. Diseases associated with PLIN3 include sebaceous adenoma, and vaccinia, and among its related super-pathways are Clathrin-dependent protein traffic. GO annotations related to this gene include protein binding. An important paralog of this gene is PLIN2.
Immunogen	A synthetic peptide corresponding to amino acids 1-100 of mouse Plin3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µl
Buffer	In Tris-glycine, 150 mM NaCl (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C or -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Plin3 perilipin 3 [Mus musculus (house mouse)]
------------------	--

Official Symbol	PLIN3
Synonyms	PLIN3; perilipin 3; Tip47; M6prbp1; 1300012C15Rik; perilipin-3; cargo selection protein TIP47; mannose-6-phosphate receptor binding protein 1; mannose-6-phosphate receptor-binding protein 1;
Entrez Gene ID	66905
Protein Refseq	NP_080112
UniProt ID	Q9DBG5
Chromosome Location	17; 17 D