



Anti-CH2 (N-terminal) polyclonal antibody (DPAB-DC2978)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus cat Fel D1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Cat
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Liquid
Size	100 µg
Buffer	In buffer containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CH2 major allergen I, polypeptide chain 2 [Felis catus (domestic cat)]
Official Symbol	CH2
Synonyms	CH2; Fel d I; major allergen I polypeptide chain 2; AG4; allergen FdI; allergen Cat-1; allergen Fel d I-B;

Entrez Gene ID	677879
Protein Refseq	NP_001041619
UniProt ID	P30440
Chromosome Location	chromosome: E2