



Rabbit anti-Mouse FBXO28 polyclonal antibody (DPAB-DC2988)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FBXO28 (F-box protein 28) is a protein-coding gene.
Specificity	Specific to recombinant protein GX0462. This antibody detects mFBXO28 protein.
Immunogen	Recombinant GST fusion protein corresponding to 104 amino acids of mouse Fbxo28.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Fbxo28 F-box protein 28 [Mus musculus (house mouse)]
-----------	--

Official Symbol	FBXO28
Synonyms	FBXO28; F-box protein 28; Fbx28; mKIAA0483; D1Ert578e; 4833428J17Rik; 5730505P19Rik; F-box only protein 28;
Entrez Gene ID	67948
Protein Refseq	NP_780336
UniProt ID	Q8BIG4
Chromosome Location	1 H5; 1 84.93 cM
Function	molecular_function;