



Anti-CRISP2 (aa 154-243) polyclonal antibody (DPAB-DC3106)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CRISP2 (cysteine-rich secretory protein 2) is a protein-coding gene. Diseases associated with CRISP2 include epididymitis, and male infertility. An important paralog of this gene is CLEC18B.
Immunogen	CRISP2 (NP_003287, 154 a.a. ~ 243 a.a) partial recombinant protein with GST tag. The sequence is YCPNQDSLKYYYVCQYCPAGNNMNRKNTPTYQQGTPCAGCPDDCDKGLCTNSCQYQDLLSN CDSLKNTAGCEHELLKEKCKATCLCENKIY
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [CRISP2 cysteine-rich secretory protein 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	CRISP2
Synonyms	CRISP2; cysteine-rich secretory protein 2; CT36; TPX1; TSP1; GAPDL5; CRISP-2; cancer/testis antigen 36; testis-specific protein TPX-1; testis specific protein 1 (probe H4-1 p3-1); glyceraldehyde-3-phosphate dehydrogenase-like 5;
Entrez Gene ID	7180
Protein Refseq	NP_001135879
UniProt ID	A0A024RD74
Chromosome Location	6p12.3