



Anti-CEV1 (N-terminal) polyclonal antibody (DPAB-DC3399)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of Arabidopsis thaliana CEV1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CEV1 cellulose synthase A catalytic subunit 3 [UDP-forming] [Arabidopsis thaliana (thale cress)]
Official Symbol	CEV1

Synonyms	CEV1; ATCESA3; ATH-B; CELLULOSE SYNTHASE 3; CESA3; CEV1; CONSTITUTIVE EXPRESSION OF VSP 1; ISOXABEN RESISTANT 1; IXR1; K2A11.4; K2A11_4; cellulose synthase A catalytic subunit 3 [UDP-forming];
Entrez Gene ID	830399
Protein Refseq	NP_196136
UniProt ID	Q941L0
Chromosome Location	chromosome: 5
Pathway	cellulose biosynthesis;
Function	cellulose synthase activity; transferase activity, transferring glycosyl groups;