



# Anti-HIST1H2BH (internal region) polyclonal antibody (DPAB-DC3414)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a member of the histone H2B family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6.
----------------------------	---

<b>Specificity</b>	This antibody detects endogenous levels of total HIST1H2BH protein.
<b>Immunogen</b>	A synthetic peptide corresponding to internal of human HIST1H2BH.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB (Cell lysate), IHC-P, IF, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
<b>Preservative</b>	0.02% Sodium Azide

**Storage**

Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

Gene Name	<a href="#">HIST1H2BH histone cluster 1, H2bh [ Homo sapiens (human) ]</a>
Official Symbol	HIST1H2BH
Synonyms	HIST1H2BH; histone cluster 1, H2bh; H2B/j; H2BFJ; histone H2B type 1-H; histone H2B.j; histone 1, H2bh; H2B histone family, member J;
Entrez Gene ID	<a href="#">8345</a>
Protein Refseq	<a href="#">NP_003515</a>
UniProt ID	<a href="#">Q93079</a>
Chromosome Location	6p21.3
Pathway	Alcoholism; Amyloids; Cell Cycle, Mitotic; Cellular responses to stress.
Function	DNA binding; protein heterodimerization activity;