



Magic™ Anti-NSRP1 (C-terminal) polyclonal antibody (DPAB-DC3460)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NSRP1 (nuclear speckle splicing regulatory protein 1) is a protein-coding gene. Diseases associated with NSRP1 include renal cell carcinoma. GO annotations related to this gene include protein binding and mRNA binding.
Specificity	At least three isoforms are known to exist; this antibody will detect all three.
Immunogen	A synthetic peptide corresponding to 16 amino acids at C-terminus of human NSRP1.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human, Mouse, Rat
Purification	Peptide affinity purification
Conjugate	Unconjugated
Applications	WB (Tissue lysate), IF, ELISA,
Format	Liquid
Concentration	1 mg/mL
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated

freezing and thawing.

GENE INFORMATION

Gene Name	NSRP1 nuclear speckle splicing regulatory protein 1 [Homo sapiens (human)]
Official Symbol	NSRP1
Synonyms	NSRP1; nuclear speckle splicing regulatory protein 1; CCDC55; NSrp70; HSPC095; coiled-coil domain containing 55; nuclear speckle-related protein 70; coiled-coil domain-containing protein 55;
Entrez Gene ID	84081
Protein Refseq	NP_001248396
UniProt ID	A0A024QZ33
Chromosome Location	17q11.2
Function	mRNA binding; poly(A) RNA binding; protein binding;