



Anti-RUNDC2A (full length) polyclonal antibody (DPAB-DC3464)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RUNDC2A (AAH41583, 1 a.a. ~ 249 a.a) full-length recombinant protein with GST tag. The sequence is MDEERSSMLPTMAAGPNSILFAINIDNKDLNGQSKFAPTVSDLLKESTQNVTLTLLKESTQG VSSVFREITASSAISILIKPEQETDPLPVVSRNVSADAKCKKERKKKKQVTNIISFDDEE DEQNSGDMFKKTPGAGESSEDNSDHSSVNIMSAFESPFGPNSNGSQSSNSWKIDSLSLNR EFGYQKLDVKSIDDEDVDENEDDVYGNSSGRKHRGHSESPEKNGAHSVTQAGVQWHDLGS LQPLPPGFR
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RUNDC2A RUN domain containing 2A [Homo sapiens (human)]
Official Symbol	RUNDC2A

Synonyms	RUNDC2A; RUN domain containing 2A; RUN domain-containing protein 2A;
Entrez Gene ID	84127
Chromosome Location	16p13.13