



Anti-NTNG2 (aa 21-119) polyclonal antibody (DPAB-DC3487)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NTNG2 (netrin G2) is a protein-coding gene. Diseases associated with NTNG2 include rett syndrome, and bipolar disorder. GO annotations related to this gene include molecular_function. An important paralog of this gene is LAMC1.
Immunogen	NTNG2 (NP_115925, 21 a.a. ~ 119 a.a) partial recombinant protein with GST tag. The sequence is ICKSWVTTDEGPTWEFYACQPKVMRLKDYVKVKVEPSGITCGDPPERFCSHENPYLCSNE CDASNPDLAHPRLMFDKEEEGLATYWQSITWSRYPSP
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [NTNG2 netrin G2 \[Homo sapiens \(human\) \]](#)

Official Symbol	NTNG2
Synonyms	NTNG2; netrin G2; Lmnt2; NTNG1; LHLL9381; bA479K20.1; netrin-G2; laminet 2; laminet-2; netrin G1; bA479K20.1 (novel protein);
Entrez Gene ID	84628
Protein Refseq	NP_115925
UniProt ID	O95898
Chromosome Location	9q34
Pathway	Axon guidance; Cell adhesion molecules (CAMs);
Function	molecular_function; protein binding;