



Anti-TRIM52 (aa 201-297) polyclonal antibody (DPAB-DC3501)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TRIM52 (tripartite motif containing 52) is a protein-coding gene. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is TRIM31.
Immunogen	TRIM52 (NP_116154, 201 a.a. ~ 297 a.a) partial recombinant protein with GST tag. The sequence is LANMVQIIRQMCPTPYRGNRSNDQGMCFKHQEALKLFCEVDKEAICVVCRESRSHKQHSV LPLEEVVQEYQEIKLETTLVGILQIEQESIHSKAYNQ
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TRIM52 tripartite motif containing 52 [Homo sapiens (human)]
Official Symbol	TRIM52

Synonyms	TRIM52; tripartite motif containing 52; RNF102; tripartite motif-containing protein 52; RING finger protein 102; tripartite motif-containing 52;
Entrez Gene ID	84851
Protein Refseq	NP_116154
UniProt ID	L0CQ38
Chromosome Location	5q35.3
Function	zinc ion binding;