



Anti-SPRYD3 (aa 1-100) polyclonal antibody (DPAB-DC3507)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SPRYD3 (SPRY domain containing 3) is a protein-coding gene. Diseases associated with SPRYD3 include agammaglobulinemia. An important paralog of this gene is RANBP10.
Immunogen	SPRYD3 (NP_116229, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. The sequence is MRRTRRPRFVLMNKMDDLNLHYRFLNWRRRIREIREVRAFRYQERFKHILVDGDTLSYHG NSGEVGCYVASRPLTKDSNYFEVSIVDSGVRGTIAVGLVP
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SPRYD3 SPRY domain containing 3 [Homo sapiens (human)]
Official Symbol	SPRYD3

Synonyms	SPRYD3; SPRY domain containing 3; SPRY domain-containing protein 3;
Entrez Gene ID	84926
Protein Refseq	NP_116229
UniProt ID	A0A024RAX4
Chromosome Location	12q13.13