



Anti-RNQ1 polyclonal antibody (DPAB-DC3514)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with <i>S. cerevisiae</i> RNQ1.
Immunogen	A synthetic peptide corresponding to <i>Saccharomyces cerevisiae</i> RNQ1. The sequence is SQQNNNGNQNRV
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RNQ1 Rnq1p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	RNQ1

Synonyms	RNQ1; Rnq1p;
Entrez Gene ID	850329
Protein Refseq	NP_009902
UniProt ID	P25367
Chromosome Location	chromosome: III
Function	molecular_function;