



Anti-TAF2 (C-terminal) polyclonal antibody (DPAB-DC3520)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TAF2 (TAF2 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 150kDa) is a protein-coding gene. Diseases associated with TAF2 include colon adenocarcinoma, and adenocarcinoma, and among its related super-pathways are Assembly of RNA Polymerase-II Initiation Complex and Apoptotic Pathways in Synovial Fibroblasts. GO annotations related to this gene include metalloproteinase activity and transcription regulatory region DNA binding.
Immunogen	A synthetic peptide corresponding to C-terminus of <i>Saccharomyces cerevisiae</i> TAF2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TAF2 Taf2p [Saccharomyces cerevisiae S288c]
Official Symbol	TAF2
Synonyms	TAF2; TAF150; TSM1; Taf2p;
Entrez Gene ID	850409
Protein Refseq	NP_009971
UniProt ID	P23255
Chromosome Location	chromosome: III
Pathway	Basal transcription factors;
Function	contributes_to RNA polymerase II core promoter sequence-specific DNA binding transcription factor activity involved in preinitiation complex assembly; chromatin binding;