



Anti-GSP1 polyclonal antibody (DPAB-DC3524)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Yeast Gsp1p (Ran).
Immunogen	Recombinant GST fusion protein corresponding to Yeast GSP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate), IF,
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	GSP1 Gsp1p [Saccharomyces cerevisiae S288c]
Official Symbol	GSP1
Synonyms	GSP1; CNR1; CST17; Gsp1p;

Entrez Gene ID	851000
Protein Refseq	NP_013396
UniProt ID	P32835
Chromosome Location	chromosome: XII
Pathway	RNA transport; Ribosome biogenesis in eukaryotes;
Function	GTP binding; GTPase activity; nucleotide binding;