



Anti-RPN5 polyclonal antibody (DPAB-DC3526)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody react with <i>S. cerevisiae</i> RPN5.
Immunogen	Recombinant protein corresponding to <i>Saccharomyces cerevisiae</i> RPN5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate), IP,
Format	Liquid
Size	100 µl
Buffer	In 10 mM Tris-HCl, 100 mM NaCl, pH 7.4 (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RPN5 proteasome regulatory particle lid subunit RPN5 [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	RPN5
Synonyms	RPN5; NAS5; proteasome regulatory particle lid subunit RPN5;

Entrez Gene ID	851408
Protein Refseq	NP_010134
UniProt ID	Q12250
Chromosome Location	chromosome: IV
Pathway	AUF1 (hnRNP D0) destabilizes mRNA; Antigen processing-Cross presentation; Proteasome degradation; CDK-mediated phosphorylation and removal of Cdc6
Function	molecular_function;