



Anti-STN1 (N-terminal) polyclonal antibody (DPAB-DC3529)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of <i>Saccharomyces cerevisiae</i> STN1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	STN1 Stn1p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	STN1
Synonyms	STN1; Stn1p;

Entrez Gene ID	851655
Protein Refseq	NP_010367
UniProt ID	P38960
Chromosome Location	chromosome: IV
Function	single-stranded telomeric DNA binding;