



Anti-BMH2 (C-terminal) polyclonal antibody (DPAB-DC3532)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to C-terminus of <i>Saccharomyces cerevisiae</i> BMH2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	BMH2 Bmh2p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	BMH2
Synonyms	BMH2; SCD3; Bmh2p;

Entrez Gene ID	851676
Protein Refseq	NP_010384
UniProt ID	P34730
Chromosome Location	chromosome: IV
Pathway	Cellular response to heat stress; HSF1 activation.
Function	DNA replication origin binding; phosphoserine binding; protein domain specific binding;