



Anti-RUB1 (full length) polyclonal antibody (DPAB-DC3535)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length <i>Saccharomyces cerevisiae</i> RUB1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Lyophilized
Size	500 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot after reconstitution to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RUB1 NEDD8 family protein RUB1 [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	RUB1

Synonyms	RUB1; NEDD8 family protein RUB1;
Entrez Gene ID	851717
Protein Refseq	NP_010423
UniProt ID	Q03919
Chromosome Location	chromosome: IV
Function	protein tag;