



Anti-SUP35 (aa 494-507) polyclonal antibody (DPAB-DC3538)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with <i>S. cerevisiae</i> SUP35.
Immunogen	A synthetic peptide corresponding to amino acids 494-507 of <i>Saccharomyces cerevisiae</i> SUP35.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Transfected lysate),
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [SUP35 Sup35p \[*Saccharomyces cerevisiae* S288c \]](#)

Official Symbol	SUP35
Synonyms	SUP35; GST1; PNM2; SAL3; SUF12; SUP2; SUP36; Sup35p;
Entrez Gene ID	851752
Protein Refseq	NP_010457
UniProt ID	P05453
Chromosome Location	chromosome: IV
Pathway	Eukaryotic Translation Termination; Metabolism of proteins; Nonsense Mediated Decay (NMD) independent of the Exon Junction Complex (EJC); Translation.
Function	GTP binding; GTPase activity; mRNA binding; nucleotide binding