



Anti-RAD9 (aa 1249-1263) polyclonal antibody (DPAB-DC3542)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No reactivity is expected against the human analog of RAD9.
Immunogen	A synthetic peptide corresponding to amino acids 1249-1263 residues surrounding S1260 of <i>Saccharomyces cerevisiae</i> RAD9.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [RAD9 Rad9p \[*Saccharomyces cerevisiae* S288c \]](#)

Official Symbol	RAD9
Synonyms	RAD9; Rad9p;
Entrez Gene ID	851803
Protein Refseq	NP_010503
UniProt ID	P14737
Chromosome Location	chromosome: IV
Pathway	Cell Cycle and Cell Division; Cell cycle - yeast; Meiosis - yeast;
Function	double-stranded DNA binding; enzyme activator activity; histone binding;