



Anti-SPT3 (aa 1-200) polyclonal antibody (DPAB-DC3544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids 1-200 <i>Saccharomyces cerevisiae</i> SPT3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	SPT3 Spt3p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	SPT3
Synonyms	SPT3; Spt3p;

Entrez Gene ID	852001
Protein Refseq	NP_010680
UniProt ID	P06844
Chromosome Location	chromosome: IV
Function	protein heterodimerization activity; transcription cofactor activity;