



Anti-RAD30 (aa 284-300) polyclonal antibody (DPAB-DC3545)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to amino acids 284-300 of <i>Saccharomyces cerevisiae</i> RAD30.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In buffer containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RAD30 Rad30p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	RAD30
Synonyms	RAD30; DBH1; Rad30p;

Entrez Gene ID	852028
Protein Refseq	NP_010707
UniProt ID	Q04049
Chromosome Location	chromosome: IV
Pathway	DNA Damage Bypass; Translesion synthesis by DNA polymerases bypassing lesion on DNA template.
Function	DNA binding; DNA-directed DNA polymerase activity; damaged DNA binding; metal ion binding