



Anti-PKH1 (aa 750-765) polyclonal antibody (DPAB-DC3546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 750-765 of yeast PKH1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PKH1 Pkh1p [Saccharomyces cerevisiae S288c]
Official Symbol	PKH1
Synonyms	PKH1; Pkh1p;

Entrez Gene ID	852101
Protein Refseq	NP_010778
UniProt ID	Q03407
Chromosome Location	chromosome: IV
Pathway	Activation of PKB; CD28 co-stimulation; Constitutive PI3K/AKT Signaling in Cancer; DAP12 interactions.
Function	ATP binding; kinase activity; nucleotide binding; protein kinase activity