



Anti-TFG2 (full length) polyclonal antibody (DPAB-DC3555)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length <i>Saccharomyces cerevisiae</i> TFG2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	TFG2 Tfg2p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	TFG2
Synonyms	TFG2; Tfg2p;

Entrez Gene ID	852888
Protein Refseq	NP_011519
UniProt ID	P41896
Chromosome Location	chromosome: VII
Pathway	Basal transcription factors; Eukaryotic Transcription Initiation; Gene Expression; RNA Pol II CTD phosphorylation and interaction with CE.
Function	ATP binding; DNA binding; catalytic activity; core RNA polymerase II binding transcription factor activity