



Anti-NAS6 polyclonal antibody (DPAB-DC3559)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody react with <i>S. cerevisiae</i> NAS6.
Immunogen	Recombinant protein corresponding to <i>Saccharomyces cerevisiae</i> NAS6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate), IP,
Format	Liquid
Size	100 µl
Buffer	In PBS (50% glycerol, 1 mg/mL BSA, 0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	NAS6 Nas6p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	NAS6
Synonyms	NAS6; Nas6p;

Entrez Gene ID	853147
Protein Refseq	NP_011748
UniProt ID	P50086
Chromosome Location	chromosome: VII
Pathway	AUF1 (hnRNP D0) destabilizes mRNA; Antigen processing-Cross presentation; Proteasome degradation; CDK-mediated phosphorylation and removal of Cdc6
Function	molecular_function;