



Anti-GCN5 (aa 1-300) polyclonal antibody (DPAB-DC3560)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids 1-300 <i>Saccharomyces cerevisiae</i> GCN5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	GCN5 Gcn5p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	GCN5
Synonyms	GCN5; AAS104; ADA4; SWI9; Gcn5p;

Entrez Gene ID	853167
Protein Refseq	NP_011768
UniProt ID	Q03330
Chromosome Location	chromosome: VII
Function	H3 histone acetyltransferase activity; N-acetyltransferase activity; histone acetyltransferase activity; lysine-acetylated histone binding