



Anti-MAD2 polyclonal antibody (DPAB-DC3562)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Yeast Mad2p.
Immunogen	Recombinant His6 fusion protein corresponding to Yeast MAD2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate), IF,
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	MAD2 Mad2p [Saccharomyces cerevisiae S288c]
Official Symbol	MAD2
Synonyms	MAD2; Mad2p;

Entrez Gene ID	853422
Protein Refseq	NP_012504
UniProt ID	P40958
Chromosome Location	chromosome: X
Pathway	Cell Cycle and Cell Division; Cell cycle - yeast; Meiosis - yeast;
Function	molecular_function;