



Anti-GAL11 (aa 1-300) polyclonal antibody (DPAB-DC3570)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids 1-300 <i>Saccharomyces cerevisiae</i> GAL11.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	GAL11 Gal11p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	GAL11
Synonyms	GAL11; ABE1; MED15; RAR3; SDS4; SPT13; Gal11p;

Entrez Gene ID	854106
Protein Refseq	NP_014591
UniProt ID	P19659
Chromosome Location	chromosome: XV
Function	RNA polymerase II activating transcription factor binding; RNA polymerase II transcription coactivator activity involved in preinitiation complex assembly; RNA polymerase II transcription cofactor activity; RNA polymerase II transcription factor recruitin