



Anti-TOA1 (full length) polyclonal antibody (DPAB-DC3573)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length <i>Saccharomyces cerevisiae</i> TOA1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	TOA1 Toa1p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	TOA1
Synonyms	TOA1; Toa1p;

Entrez Gene ID	854369
Protein Refseq	NP_014837
UniProt ID	P32773
Chromosome Location	chromosome: XV
Pathway	Basal transcription factors; Gene Expression; RNA Polymerase II Promoter Escape; RNA Polymerase II Transcription Initiation.
Function	contributes_to TBP-class protein binding; contributes_to TBP-class protein binding RNA polymerase II transcription factor activity involved in preinitiation complex assembly;