



Anti-YPK2 (N-terminal) polyclonal antibody (DPAB-DC3578)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of yeast YPK2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	YPK2 Ypk2p [Saccharomyces cerevisiae S288c]
Official Symbol	YPK2
Synonyms	YPK2; YKR2; Ypk2p;

Entrez Gene ID	855130
Protein Refseq	NP_013822
UniProt ID	P18961
Chromosome Location	chromosome: XIII
Function	ATP binding; kinase activity; nucleotide binding; protein kinase activity