



Anti-RFA2 (aa 115-129) polyclonal antibody (DPAB-DC3584)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 115-129 of <i>Saccharomyces cerevisiae</i> RFA2. The sequence is SSGKGYGSQVAQQFE
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate), IP, ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RFA2 Rfa2p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	RFA2

Synonyms	RFA2; BUF1; RPA2; RPA32; Rfa2p;
Entrez Gene ID	855404
Protein Refseq	NP_014087
UniProt ID	P26754
Chromosome Location	chromosome: XIV
Pathway	DNA Repair; DNA replication; Dual incision reaction in GG-NER; Global Genomic NER (GG-NER)
Function	double-stranded DNA binding; nucleic acid binding; sequence-specific DNA binding; single-stranded DNA binding