



Anti-ULP1 polyclonal antibody (DPAB-DC3586)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to <i>Saccharomyces cerevisiae</i> ULP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IF, ELISA,
Format	Lyophilized
Size	500 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot after reconstitution to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ULP1 Ulp1p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	ULP1
Synonyms	ULP1; NIB1; Ulp1p;

Entrez Gene ID	856087
Protein Refseq	NP_015305
UniProt ID	Q02724
Chromosome Location	chromosome: XVI
Function	SUMO-specific protease activity; cysteine-type peptidase activity; hydrolase activity; peptidase activity