



Anti-SUA7 (full length) polyclonal antibody (DPAB-DC3587)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length <i>Saccharomyces cerevisiae</i> SUA7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	SUA7 Sua7p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	SUA7
Synonyms	SUA7; SOH4; Sua7p;

Entrez Gene ID	856201
Protein Refseq	NP_015411
UniProt ID	P29055
Chromosome Location	chromosome: XVI
Pathway	Basal transcription factors; Eukaryotic Transcription Initiation; RNA Polymerase II Pre-transcription Events; RNA Polymerase II Transcription.
Function	RNA polymerase II core binding; TBP-class protein binding; TBP-class protein binding RNA polymerase II transcription factor activity involved in preinitiation complex assembly; core RNA polymerase II recruiting transcription factor activity