



Anti-SCH9 (N-terminal) polyclonal antibody (DPAB-DC3590)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of yeast SCH9.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SCH9 Sch9p [Saccharomyces cerevisiae S288c]
Official Symbol	SCH9
Synonyms	SCH9; HRM2; KOM1; Sch9p;

Entrez Gene ID	856612
Protein Refseq	NP_012075
UniProt ID	P11792
Chromosome Location	chromosome: VIII
Function	ATP binding; kinase activity; nucleotide binding; protein kinase activity