



Anti-RPN3 polyclonal antibody (DPAB-DC3591)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody react with <i>S. cerevisiae</i> RPN3.
Immunogen	Recombinant protein corresponding to <i>Saccharomyces cerevisiae</i> RPN3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate), IP,
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	RPN3 proteasome regulatory particle lid subunit RPN3 [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	RPN3
Synonyms	RPN3; SUN2; proteasome regulatory particle lid subunit RPN3;

Entrez Gene ID	856742
Protein Refseq	NP_010938
UniProt ID	P40016
Chromosome Location	chromosome: V
Pathway	AUF1 (hnRNP D0) destabilizes mRNA; Antigen processing-Cross presentation; Proteasome degradation; CDK-mediated phosphorylation and removal of Cdc6
Function	enzyme regulator activity; molecular_function;