



Anti-SRB4 (full length) polyclonal antibody (DPAB-DC3592)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length <i>Saccharomyces cerevisiae</i> SRB4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	SRB4 Srb4p [<i>Saccharomyces cerevisiae</i> S288c]
Official Symbol	SRB4
Synonyms	SRB4; MED17; Srb4p;

Entrez Gene ID	856743
Protein Refseq	NP_010939
UniProt ID	P32569
Chromosome Location	chromosome: V
Function	RNA polymerase II transcription cofactor activity; activating transcription factor binding; core RNA polymerase II recruiting transcription factor activity; structural molecule activity