



Anti-Acetyl-Lysine polyclonal antibody (DPAB-DC3948)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Detects proteins containing acetylated lysine residues in SDS-PAGE immunoblots.
Immunogen	Acetylated Lysine conjugated with KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB (Tissue lysate), IF, IP, ELISA,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.0 (50% glycerol)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Acetyllysine (or acetylated lysine) is an acetyl-derivative of the amino acid lysine. There are multiple forms of acetyllysine - this article refers to N-ε-acetyl-L-lysine. The other form is N-α-acetyl-L-lysine. In proteins, the acetylation of lysine resi
---------------------	---

Keywords

Acetyllysine; Acetylated Lysine
