



Anti-Activin A polyclonal antibody (DPAB-DC3954)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to human activin A. The sequence is GLECDGKVNICKKKQFFVSFKDIGWNDWIIAPSGYHANYCEGECPSHIAGTSGSSLSFHS TVINHYRMRGHSPFANLKSCCVPTKLRPMSMLYYDDGQNIKKDIQNMIVEECGCS
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purification
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from PBS
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with 100 uL distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords Active A