



Anti-Alanine-beta polyclonal antibody (DPAB-DC3962)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	beta-alanine conjugated with carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ICC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with deionized water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	β-Alanine (or beta-alanine) is a naturally occurring beta amino acid, which is an amino acid in which the amino group is at the β-position from the carboxylate group (i.e., two atoms away, see Figure 1). The IUPAC name for β-alanine is 3-aminopropanoic ac
---------------------	---

KeywordsAlanine-beta
