



Anti-Aldehyde Oxidase (native) polyclonal antibody (DPAB-DC3964)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified rat aldehyde oxidase.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB (Tissue lysate),
Format	Liquid
Size	25 µg
Buffer	In PBS (0.1% proclin, 2.0% Block Ace)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Aldehyde oxidase (AO) is a metabolizing enzyme, located in the cytosolic compartment of tissues in many organisms. AO catalyzes the oxidation of aldehydes into carboxylic acid, and in addition, catalyzes the hydroxylation of some heterocycles. It can also
Keywords	AO; Aldehyde Oxidase