



Anti-Aldosterone (C-terminal) polyclonal antibody (DPAB-DC3965)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	4-PREGNEN-11b, 21-DIOL-3, 18, 20-TRIONE 3-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Aldosterone is a steroid hormone (mineralocorticoid family) produced by the outer section (zona glomerulosa) of the adrenal cortex in the adrenal gland. It plays a central role in the regulation of blood pressure mainly by acting on the distal tubules and
Keywords	Aldosterone